



Recombinant Shikimate dehydrogenase (aroE)

Product Code	CSB-MP719818SMZ
Abbreviation	aroE
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XAV7
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MSERLSGHTL LVSLLATPIR HSLSPKMNE AYAKLGLDYA YLAFEVGT EQ LADAVQGIRA LGIRGSNVSM PNKEAILPLL DDLSPAAELV GAVNTVVNKD GKGHLVGHIT DGIGALRALA DEGVS VKNKI IFLAGVGGAG KAIQVQLAFD GAKEIRLFNR QATRLSSVQK LVTKLNQLTR TKVTLQDLED QTAFKEAIRE SHLFIDATSV GMKPLENLSL ITDPELIRPD LVVFDIVYSP AETKLLAFAR QHGAQKQVING LGMVLYQGAE AFKLITGQDM PVDAIKPLLG DK</pre>
Source	Mammalian cell
Target Names	aroE
Protein Names	Recommended name: Shikimate dehydrogenase EC= 1.1.1.25
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>