



Recombinant Uridylate kinase (pyrH)

Product Code	CSB-BP719835SMZ
Abbreviation	pyrH
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q5XDH4
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MEPKYQRILI KLSGEALAGE KGVGIDIPTV QAIAKEIAEV HVSGVQIALV IGGGNLWRGE PAADAGMDRV QADYTGMLGT VMNALVMADS LQHYGVDTRV QTAIPMQNVA EPYIRGRALR HLEKNRIVVF GAGIGSPYFS TDTTAALRAA EIEADAILMA KNGVDGVYNA DPKKDANAVK FDELTHGEVI KRGLKIMDAT ASTLSMDNDI DLVVFNMNEA GNIQRVVFGE HIGTTVSNKV CD
Source	Baculovirus
Target Names	pyrH
Protein Names	<p>Recommended name: Uridylate kinase Short name= UK EC= 2.7.4.22</p> <p>Alternative name(s): Uridine monophosphate kinase Short name= UMP kinase</p> <p>Short name= UMPK</p>
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>