



Recombinant Segregation and condensation protein A (scpA)

Product Code	CSB-EP719837SMZ
Abbreviation	scpA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XDP4
Product Type	Recombinant Protein
Immunogen Species	Streptococcus pyogenes serotype M6 (strain ATCC BAA-946 / MGAS10394)
Purity	≥85% (SDS-PAGE)
Sequence	MDIKLKDFEG PLDLLLHLVS QYKVDIYEV IVEVIEQYLN YIETLQVMKL EVAGDYMLMA SQLMLIKSRR LLPKVVEHIE EEDLEQDLLE KIEEYSRFKA VSQALAKQHD QRAKWYSKPK QELIFEDAIL QEDKTVMDLF LAFSNIMAAK RAVLKNNHTV IERDDYKIED MMASIKQRLE KESVISLSAI FEECQTLNEV ISIFLASLEL IKLHVVFVEQ LSNFGAILR KEKK
Source	E.coli
Target Names	scpA
Protein Names	Recommended name: Segregation and condensation protein A
Expression Region	1-234
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.