



Recombinant Rat Probable tRNA pseudouridine synthase 2 (Trub2)

Product Code	CSB-YP719852RA
Abbreviation	Trub2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5XFW2
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MGSSGLARLQ GLFAVYKPPG LKWLHVRETV ELQLLKGLNA QQPSAPEQHV RFLGPEVEGS EEKTLTLTAT SVPSLTTHRL VRGPAFRNLK IGVGHRLDVQ ASGVLVLAVG HGRSLLTDMY DAHLTKDYTV RGLLGKATDN FCEDGQLIEK TTYDHSVTRER LDRILAMIQG SHQKALVMYS NLDLKSQEAY ERAVQGVIRP MNKSPMLIAG IRCLHFAPPE FLLEVQCMHE TQQQLRRLVH EIGLELKTTA VCMQVRRTRD GFFGLDDALL RTQWDLHSIQ DAIQTAAPRV ARELRKNLSL MSGHQQLPS AGQPWASRVQ APL
Source	Yeast
Target Names	Trub2
Protein Names	Recommended name: Probable tRNA pseudouridine synthase 2 EC= 5.4.99.-
Expression Region	1-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.