



# Recombinant Rat TRAF3-interacting protein 1 (Traf3ip1)

<b>Product Code</b>	CSB-MP719882RA
<b>Abbreviation</b>	Traf3ip1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5XIN3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MNAAVVRRRTQ EALGKVIRRP PLTEKLLNKP PFRYLHDIIT EVIRITGFMK GLYTDAEMKS DNVKDKDAKI SFLQKAIDVV MMVSGEPLAA KPARIVAGHE PERTNELLQL IGCCLSKLS SDEAVKRVLA GEKGDSRGRV LRSTKAQEPD NKSVKEEESR TQKEEKRSSE VKERSSSAEH KQKEELKEDS KPREKERDKE KAKEADRRH REPDRDRNRD GEREKARARA KQDRDRNND RDRETERDRD RDRRSDDGGKE KERQKDRDRD KGKDRERRKS KNGEHTRDPD REKSRDADKS EKKADISVGA SKSSTLKASK RRSKNSLEGR KEDNISAKIL DSIVSGLNDE PDQETTAPEI DDNSASLWRE NAEPEPAVKQ KGDSPSDAEG EAVPTSQDKL EVTENAEVSN ELPSSLRRIP RPGSARPAPP RVKRQUESTET LAGDRSGSGK TVSTVIIDSQ NSDNEDDEQF VVEAAPQLSE IAEIMVPSG DLEDEEKHGG LVKKILETKK DYEKLQQSPK PGEKERSLIF ESAWKKEKDI VSKEIEKLRV SIQTLCKSAL PLGKIMDYIQ EDVDAMQNEL QLWHSNRQH AEALSKEQSI TDSAVEPLKA ELSELEQQIK DQQDKICAVK ANILKNEEKI QKMOVHSINLS SRR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Traf3ip1
<b>Protein Names</b>	Recommended name: TRAF3-interacting protein 1
<b>Expression Region</b>	1-653
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.