



Recombinant Nocardia farcinica Succinyl-CoA ligase [ADP-forming] subunit beta (sucC)

Product Code	CSB-EP719927NAAA-B
Abbreviation	sucC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5YPP8
Product Type	Recombinant Protein
Immunogen Species	Nocardia farcinica (strain IFM 10152)
Purity	≥85% (SDS-PAGE)
Sequence	MDLFEYQAKE LFKVHGVPSS EGRVTDTAED ARAIAAEIGK PVMIKSQVKV GGRGKAGGVK YAATPEDAFT HAQNILGLDI KGHITKKILV AEAKDIAEEY YISFLLDRAN RTYLAMCSVE GGMEIEEVAA TKPERLAKVP VDAVKGVDLA FARSIAEQGH LPAEVLDAAA VTIQKLWEVF VAEDATLVEV NPLVRTPQDE ILALDGKVTI DENAGFRHPE HDEFADVDAT DPLELKAKEN DLNYVKLDGE VGIIGNGAGL VMSTLDVVAY AGENHGGVKP ANFLDIGGGA SAEVMAAGLD VILNDAQVKS VFNVFEGGIT ACDVANGIV KALEILGDEA NKPLVVRLDG NKVDEGRKIL VEAHPLVTL AQTMDGADK AAELAAAR
Source	E.coli
Target Names	sucC
Protein Names	Recommended name: Succinyl-CoA ligase [ADP-forming] subunit beta EC=6.2.1.5 Alternative name(s): Succinyl-CoA synthetase subunit beta Short name=SCS-beta
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.