



Recombinant *Oryza sativa* subsp. japonica Probable calcium-binding protein CML29 (CML29)

Product Code	CSB-YP719970OFG
Abbreviation	CML29
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5Z676
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MATPAGGRPS QPQPQAQQLS VDFEALSYIS SLVEAFQAFD SDNDGLVTAP ELRGLLASLG LDKPEHEVRD MLARADADRD GKLSVEELLD VMNAGQLGLG ALGALLQSAV PALESAAGPD GVLGADELAR LLSVMGTASV EDCMEIACM DGDGDGAISV EEFRLMAQLL
Source	Yeast
Target Names	CML29
Protein Names	Recommended name: Probable calcium-binding protein CML29 Alternative name(s): Calmodulin-like protein 29
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.