



Recombinant Chicken SPRY domain-containing protein 7 (SPRYD7)

Product Code	CSB-MP719981CH
Abbreviation	SPRYD7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZHV7
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MAASVFCCLR WCRDGGAGHI PLKEMPAVQL DTQRMGTDVV IVKNGRRICG TGGCLANAPL HQNKSYFEFK IQSTGIWGIG VATQKANLNQ IPLGRDVHSL VMRNDGALYY NNEEKNRLPA NSLPQEGDVV GITYDHVELN VYLNKGNMHC PASGIRGTVY PVVYVDDSAI LDCQFSEFYH TPPPGEKIL FEQQIF
Source	Mammalian cell
Target Names	SPRYD7
Protein Names	Recommended name: SPRY domain-containing protein 7 Alternative name(s): Chronic lymphocytic leukemia deletion region gene 6 protein homolog Short name= CLL deletion region gene 6 protein homolog
Expression Region	1-196
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.