



Recombinant Chicken Amyloid beta A4 precursor protein-binding family B member 1-interacting protein (APBB1IP)

Product Code	CSB-MP720019CH
Abbreviation	APBB1IP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZL23
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	MEQTCDDIDE MFSNLLGEMD MLTQSLGVET VQTPSPKVTN EEFSFTVGFK DLNESLNALE DKDLDALMAD LVADINEVEQ RTLQAQKTSG NQQSVVTQPS TGTNNDFC SK LSPCATITGQ FKNDLPPPPP APDLDLPPPP PPPPEPLSQ EEQEARAKAD KIKLALEK LK EAKIKKLVVK VHMYDNSTKS LMVDERQVTR DVL DNLFEKT HCDCSVDWCL YEVYPELQIE RFFEDHENVV EVLSDWTRDS ENKVL FLEKK EKYALFKNPQ NYLANKGKN ESKE MNDKSK EALLEESFCG ASVIVPELEG ALY LKEDGKK SWKRRYFLLR ASGIYYVPKG KTKTSRDLMC FIQFENMNVY YGSQHKVKYK APTDHC FVLK HPQIQKESQY IKYLCCDDRA TLHQWVTGIR IAKYGK TLYD NYKCAVKKAG LSSQWANQGT LEPAAPTGSL SAGAVQANGQ IPRVVL PSSA EVAETQKKVD PAIRA
Source	Mammalian cell
Target Names	APBB1IP
Protein Names	Recommended name: Amyloid beta A4 precursor protein-binding family B member 1-interacting protein Alternative name(s): APBB1-interacting protein 1
Expression Region	1-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.