



Recombinant Chicken Methyltransferase-like protein 9 (METTL9)

Product Code	CSB-BP720041CH
Abbreviation	METTL9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZMH6
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	>85% (SDS-PAGE)
Sequence	GPARYL RSPLSRSLYA NMMGSHGPPA PGAGENHQWY VCNTEQLSES LQPIFVQSYL DQGTQIFLNN SIEKSGWLFQ QLYHSFVSSI FSLFMSRTSI NGLLGRGSMF VFSPEQFQRL LKINPEWKSH RLLDLGAGDG EVTKVMSPHF EEIYATELSE TMIWQLQKKK YRVLGINEWQ NTGFQYDVIS CLNLLDRCDQ PLTVLKDTRS VLEPTRGRVI LALVLPFHPY VENVGGKWEK PSEVLEIKGH TWEEQVNSLP EVFGKAGFAI EAFTRLPYLC EGDMYNDYYV LDDAVFVLKP V
Source	Baculovirus
Target Names	METTL9
Protein Names	Recommended name: Methyltransferase-like protein 9
Expression Region	25-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.