



Recombinant Legionella pneumophila subsp. pneumophila Deoxycytidine triphosphate deaminase (dcd)

Product Code	CSB-YP720046LDL
Abbreviation	dcd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZRC6
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	>85% (SDS-PAGE)
Sequence	MSIKSDRWIE KMALEHGMIS PFQAGQVREN QNGRIISYGV SSYGYDVRCS NEFKIFTNIN SAIVDPKAFD ENSFVDVQSD VCIIPPNSFA LARTVEYFRI PRNILTICLG KSTYARCGII VNVTPLEPEW EGHVTLEFSN TTTLPKIYA YEGVAQMLFL EANEVCAVS YRDRNGKYQQG TGVTLPR T
Source	Yeast
Target Names	dcd
Protein Names	Recommended name: Deoxycytidine triphosphate deaminase Short name= dCTP deaminase EC= 3.5.4.13
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.