



Recombinant Legionella pneumophila subsp. pneumophila Ribosomal RNA small subunit methyltransferase G (rsmG)

Product Code	CSB-EP720047LDL
Abbreviation	rsmG
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5ZRJ1
Product Type	Recombinant Protein
Immunogen Species	Legionella pneumophila subsp. pneumophila (strain Philadelphia 1 / ATCC 33152 / DSM 7513)
Purity	≥85% (SDS-PAGE)
Sequence	MIDNTKIRSL LEKGLAEFKL DSIGDLLLDF LLLLNKWNKT YNLTAIRDIE TMVSKHLFDS LAILPWIKGN HIIDVGTGPG LPGIPLAIK PDLQFVLLDS NGKKIRFLNE VKRQLNIKNI EPIQIRVENY HPNQGFDTVI SRAFSSLEQM IKWTQHLVAQ DGLWLAMKGR FPDTELVPIH QTYRVERYAV PGIEGERCCV LINNTNKE
Source	E.coli
Target Names	rsmG
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase G EC= 2.1.1.170 Alternative name(s): 16S rRNA 7-methylguanosine methyltransferase Short name= 16S rRNA m7G methyltransferase
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.