



Recombinant Methanosarcina mazei V-type ATP synthase subunit E (atpE)

Product Code	CSB-BP720101MSN
Abbreviation	atpE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60183
Product Type	Recombinant Protein
Immunogen Species	Methanosarcina mazei (strain ATCC BAA-159 / DSM 3647 / Goe1 / Go1 / JCM 11833 / OCM 88) (Methanosarcina frisia)
Purity	>85% (SDS-PAGE)
Sequence	MGLEIVVKDI QEGARAEVSR IKAEGDAKAS EIINEAKEIQ KKTGLGDSLAK AEEDLQSLHQ QVISSANLEV KRITLNKRKE LLDTVYNQTV ENIKSMPASK KEELLKSILD KHEASGARAY SSKESEELVK KLTSLSYAGN IDSIGGIVLE NEDRTVRLDF TYDSILKSVY ERSLKQISDI LYG
Source	Baculovirus
Target Names	atpE
Protein Names	Recommended name: V-type ATP synthase subunit E Alternative name(s): V-ATPase subunit E
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.