



Recombinant Methylococcus capsulatus tRNA (cmo5U34)-methyltransferase (cmoA)

Product Code	CSB-MP720156MFM
Abbreviation	cmoA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q605H3
Product Type	Recombinant Protein
Immunogen Species	Methylococcus capsulatus (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MSRDEIYRQS QSFLDDFGFG ESVATVFDDM LERSVPFYAE LQRMIAEMAG DFAMPHTRIY DFGCSTGTTL IGLDQAIGPR GLTLVGVDNS QEMLAKCRAK MAGHAFANAV ELIRADLNQG VVMENASLAL MILTLQFVRP LYRDRLVRAI HDGLEENGAL VLVEKVLGES SLFNRSFIKY YYEFKKRNGY TELEIAQKRE ALENLVLPYK LAENLEMLLES AGFRYVDVFF KWYNFVGIVA VK
Source	Mammalian cell
Target Names	cmoA
Protein Names	Recommended name: tRNA (cmo5U34)-methyltransferase EC= 2.1.1.-
Expression Region	1-242
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.