



Recombinant *Methylococcus capsulatus* Ribosomal RNA large subunit methyltransferase H (rlmH)

Product Code	CSB-MP720209MFM
Abbreviation	rlmH
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60BT8
Product Type	Recombinant Protein
Immunogen Species	<i>Methylococcus capsulatus</i> (strain ATCC 33009 / NCIMB 11132 / Bath)
Purity	≥85% (SDS-PAGE)
Sequence	MQLNLVCVGE RMPSWVNEGY EYARRLPRE CALMLREIAP GKRTRNCDIQ RILADEGQRM LAAVGGNAHI VTLDVGGKDW STRELAEVLA RWLREGRDIA LLVGGPDGLS DACRQRAAES WLSLRMTFPH PLVRVIVAEQ LYRAWTLLNN HPYHR
Source	Mammalian cell
Target Names	rlmH
Protein Names	Recommended name: Ribosomal RNA large subunit methyltransferase H EC=2.1.1.177 Alternative name(s): 23S rRNA (pseudouridine1915-N3)-methyltransferase 23S rRNA m3Psi1915 methyltransferase rRNA (pseudouridine-N3)-methyltransferase R
Expression Region	1-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.