



Recombinant Defective in cullin neddylation protein 1 (dcn-1)

Product Code	CSB-MP720232CXX
Abbreviation	dcn-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60YT5
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MNRLKADQKT KVGKRKQDPP ETLSEAEYDV KKCRAQNNAN LVGVSNMTSS QTVMPITHLS YASSSSQGLA QSSILSKNCP KNYENILRQF VQWTQATEPV SLNFLAKANW NIEYAMTLYY DNPNFLSSSA PATVDQSKTI QLFTQYVDKR DGLGERIGPH GMQRLIDLG YEPIDRRVLI LAWVFKAETQ CEFSLQEFTN GMASLQVDSI QGLKQKIDAL DAGMKADLTG TRDLCIFTFN YGKSAASRSL DLEMAICYWD VIFGARKPLM SQWIDFLYGQ ERMAYARLEE ELGPTNAKKV KSVWITRDTW NLFWDFIVLG KEDLSYDEE GAWPVLIDQF VDHCRENLNY PKPLRENSNA MEHQSY
Source	Mammalian cell
Target Names	dcn-1
Protein Names	Recommended name: Defective in cullin neddylation protein 1
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.