



# Recombinant Mouse Chitinase-3-like protein 1 (Chi3l1)

<b>Product Code</b>	CSB-BP720246MO
<b>Abbreviation</b>	Chi3l1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q61362
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	YKLVCFYFTS WSQYREGVGS FLPDAIQPFL CTHIYSFAN ISSDNMLSTW EWNDESNYDK LNKLKTRNTN LKTLLSVGGW KFGEKRFSEI ASNTERRTAF VRSVAPFLRS YGFDGLDLAW LYPRLRDQY FSTLIKELNA EFTKEVQPGR EKLLLSAALS AGKVAIDTGY DIAQIAQHLD FINLMTYDFH GVWRQITGHH SPLFQGQKDT RFDYRSNVNY AVQYMIRLGA QASKLLMGIP TFGKSFTLAS SENQLGAPIS GEGLPGRFTK EAGTLAYYEI CDFLKGAEVH RLSNEKVPFA TKGNQWVGYE DKESVKNKVG FLKEKLAGA MVWALDLDDF QGTCQPKEFF PLTNAIKDAL A
<b>Source</b>	Baculovirus
<b>Target Names</b>	Chi3l1
<b>Protein Names</b>	Recommended name: Chitinase-3-like protein 1 Alternative name(s): BRP39 protein Cartilage glycoprotein 39 Short name= CGP-39 Short name= GP-39
<b>Expression Region</b>	30-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.