



Recombinant Probable V-type proton ATPase subunit H 2 (vha-15)

Product Code	CSB-YP720269CXX
Abbreviation	vha-15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q619W9
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis briggsae
Purity	≥85% (SDS-PAGE)
Sequence	MAEVAHHNIP AVDMINATSR LQLEAAEVRS NKPWGSYFR SKMIQEDDYT LVTSYENAKS KDERDQVIAA NDANGQLAKT MANLITQVAK DQNVRYVLT FDDMLQEDKS RVEIFHKAH RQKRTAFS QY LGILQRQDNF IVNQMSIIA KLACFGVTRM EGQDLQYYFS FLKEQLKNST TNEYMNTTAR CLQMMLRHDE YRHEFVVDG VQTLVTALNG KTNFQLQYQL IFSVWCLTFN ADIAKKTPSL GVITALGDIL SESTKEKVIR IILASFVNIL NKVEEREIKR EAALQMVQCK TLKTLELMDA KKYDDPDLED DVKFLTEELT LSVHDLSSYD EYYSEVRSGR LTWSPVHKSE KFWRENAAKF NDKQYEVVKI LIKLLESSSD PLILCVASHD IGEYVRHYPR GKTVEQQQG KAAVMRLTA EDPNVRYHAL LAVQKLMVHN WEYLGKQLDS DAQGDGVA AK
Source	Yeast
Target Names	vha-15
Protein Names	Recommended name: Probable V-type proton ATPase subunit H 2 Short name= V-ATPase subunit H 2 Alternative name(s): Vacuolar proton pump subunit H 2
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.