



Recombinant O-phosphoseryl-tRNA (Sec) selenium transferase (secs-1)

| | |
|--------------------------|---|
| Product Code | CSB-MP720276CXX |
| Abbreviation | secs-1 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q61JN8 |
| Product Type | Recombinant Protein |
| Immunogen Species | Caenorhabditis briggsae |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKSSFGKKEG EYSRLVSKSS NKLLNSLWEK KQIPEEGWTE HTLDLFLSWL SSHDTNNRVD MIPVGAGERE GRVLTPLVQR LHSNLTHGIG RSGNLLIQQ KALGSSMLAC LSNEFAKHAL HLLGLQTVKS CIVVPLCTGM SLSLCMTSWR RRRPKAKYVI WLRIDQKSSL KSIYHAGFEA IIVEPTRDHD ALVTDVETVN RIVEQRGEEL LCVMTTSCF APRSPDNIEA ISAICAAHDV PHLVNNAYGL QSEETIRKIA AAHECGRVDA VVQSLDKNFQ VPVGGALIAG FKQSHIQSIA QAYPGRASSV PSRDLVLTFL YQQQSAFLEP FQKQKQMFLK MRRKLTSFAE NVGECVYDVP ENEISLAMTL STIPPTKQTL FGSVLFSGRI TGARVVQSSQ SKTTIEGCEF VNFGSHTAEQ HGGYLNIACS IGMADHELEE LFTRLTSSYA KFIRQLAKED DRRGGSSGRR VPMNESFDME ND |
| Source | Mammalian cell |
| Target Names | secs-1 |
| Protein Names | Recommended name: O-phosphoseryl-tRNA(Sec) selenium transferase EC=2.9.1.2 Alternative name(s): Selenocysteine synthase Short name= Sec synthase Selenocysteinyl-tRNA(Sec) synthase Sep-tRNA:Sec-tRNA synthase Short name= |
| Expression Region | 1-482 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.