



Recombinant Rat Non-muscle caldesmon (Cald1)

Product Code	CSB-BP720315RA
Abbreviation	Cald1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62736
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MLSRSGSQGR RCLATLSQIA YQRNDDDEEE AARERRRRAR QERLRQKQEE ESLGQVTDQV EAHVQNSAPD EESKPATANA QVEGDDEEAAL LERLARREER RQKRLQEALE RQKEFDPTIT DGSLSVPSRR MQNNSAENET AEGEEKGESR SGRYEMEETE VVITSYQKNS YQDAEDKKKE EKEEEEEEEK LKGGNLGENQ IKDEKIKKDK EPKEEVKNFL DRKKGFTVK AQNGEFMTHK LKQTENAFSP SRSGGRASGD KEAEGAPQVE AGKRLEELRR RRGETESEEF EKLKQKQQEA ALELEELKKK REERRKVL EEQRRKQEEA DRKAREEEEEK RRLKEEIERR RAEAAEKRQK MPEDGLSEDK KPFKCFPKG SSLKIEERAE FLNKSVMQKSG VKSTHQAQAVV SKIDSRLEQY TNAIEGTKAS KPMKPAASDL PVPAEGVRNI KSMWEKGSVF SPSASGTPN KETAGLKVGV SSRINEWLTK SPDGNKSPAP KPSDLRPGDV SGKRNLWEKQ SVDKVTSPTK V
Source	Baculovirus
Target Names	Cald1
Protein Names	Recommended name: Non-muscle caldesmon Short name= CDM Alternative name(s): L-caldesmon
Expression Region	1-531
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a calmodulin- and actin-binding protein that plays an essential role in the regulation of smooth muscle and nonmuscle contraction. The conserved domain of this protein possesses the binding activities to Ca(2+)-calmodulin, actin, tropomyosin, myosin, and phospholipids. This protein is a potent inhibitor of the actin-tropomyosin activated myosin MgATPase, and serves as a mediating factor for Ca(2+)-dependent inhibition of smooth muscle contraction. Alternative splicing of this gene results in multiple transcript variants encoding distinct isoforms.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.