



# Recombinant Rat Rab GTPase-binding effector protein 2 (Rabep2)

<b>Product Code</b>	CSB-EP720321RA
<b>Abbreviation</b>	Rabep2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q62835
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAAAPAALAL DPQPQEEQKD ASESSELSRL RAELAGALAE METMKAVAEV SESTKAEAVA AVQRQCQEEV ASLQAILKDS ISSYETQIAA LKQERQQQQQ DSEEKERELG HLKQLLARAH PLDSLEKQME KAHEDSEKLR EIVLPMEQEI AELKVKLLRA EELIQEIQRR PRQPASLHGS TELLPLSRNP SPPLPSEEP SGDAGPAAEA FAHNCDDSAS ISSFSLGGAA GSVSLRGPQG LSPEQEETAS LVSTGTLVPE GIYLPPPGYQ LVPDSQWEQL QVEGRQLQKE LESVRRERDE LQEGLSRNE DCAKQMQVLL AQVQNSEQLL RTLQGTVSQA QERVQRQMAE LATSHKCLSQ EVKRLNEENR GLRAEQLPSS ALQGSEQQED QDEALPSSIQ ELHQLVRHTR QQARARQQAQ EHEAERLRIE IVKLREALDE ETAAKASLEG QLRVQRETD VLEASLCSLR IETERVQQEH HKAQLTDLLS EQRKALRLQ AELETSEQVQ RDFVRLSQAL QVRLEQIRQA DTLEQVRSIL DEAPLRDIKD IKDS
<b>Source</b>	E.coli
<b>Target Names</b>	Rabep2
<b>Protein Names</b>	Recommended name: Rab GTPase-binding effector protein 2 Alternative name(s): MP13 protein Rabaptin-5beta
<b>Expression Region</b>	1-554
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.