



Recombinant Burkholderia mallei Diaminopimelate epimerase (dapF)

Product Code	CSB-MP720337BNS
Abbreviation	dapF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q62EZ3
Product Type	Recombinant Protein
Immunogen Species	Burkholderia mallei (strain ATCC 23344)
Purity	≥85% (SDS-PAGE)
Sequence	MKLSFTKMHG AGNDFVVLDTG YTRALPPLTG AQVRALADRH FGIGADQLLL VEKPTVDGAD FKYRIFNCDG GEVEHCGNGA RCFVKFVRDH GLTGKASVRV EVKHGVITLT MQDNGEVVVD MGAPVFEPAR VPFASGLDGL RREGADTLWP LPVNGVTRWI SVVSMGNPHA VQIVDDAEAF AVRVDGPAIE RDPFRPQVRN AGFMQIVSRH EVNLRVYERG AGETLACGTG ACAAVAAGIR RGRLDSPVTV HTHGGTLTIS WNGACDERAP LMMAGPATTV FEGVIELPA
Source	Mammalian cell
Target Names	dapF
Protein Names	Recommended name: Diaminopimelate epimerase Short name= DAP epimerase EC= 5.1.1.7
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.