



Recombinant Bacillus cereus Malate dehydrogenase (mdh)

Product Code	CSB-MP720507BAAD
Abbreviation	mdh
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q633K5
Product Type	Recombinant Protein
Immunogen Species	Bacillus cereus (strain ZK / E33L)
Purity	≥85% (SDS-PAGE)
Sequence	MTIKRKKVSV IGAGFTGATT AFLAQKELA DVVLVDIPQL ENPTKGKALD MLEASPVQGF DANIIGTSDY ADTADSDVVV ITAGIARKPG MSRDDLVATN SKIMKSITRD IAKHSPNAII VLTNPVDAM TYSVFKEAGF PKERVIGQSG VLDTARFRTF IAQELNFSVK DITGFVLGGH GDDMVPLVRY SYAGGIPT LIPKERLEAI VERTRKGGGE IVLLGNGSA YYAPAASLVE MTEAILKDQR RVLPAIAYLE GEYGYSPLYL GVPVILGGNG IEKIELELL ADEKEALDRS VESVRNVMKV LV
Source	Mammalian cell
Target Names	mdh
Protein Names	Recommended name: Malate dehydrogenase EC= 1.1.1.37
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.