



# Recombinant Rat C4b-binding protein beta chain (C4bpb)

<b>Product Code</b>	CSB-EP720522RA
<b>Abbreviation</b>	C4bpb
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q63515
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	LDGSC SEPPPVNNSV FVGKETEEQI LGIYLCIKGY HLVGKKSLVF DPSKEWNSTL PECLLGHCPD PVLENGKINS SGPVNISGKI MFECNDGYIL KGSNWSQCLE DHTWAPPLPI CRSRDCEPPE TPVHG YFEGE TFTSGSVVTY YCEDGYHLVG TQKVQCSDGE WSPSYPTCES IQEPPKSAEQ SALEKAILAF QESKDL CNAT ENFVRQLREG GITMEELKCS LEMKTKLKS DILLNYHS
<b>Source</b>	E.coli
<b>Target Names</b>	C4bpb
<b>Protein Names</b>	Recommended name: C4b-binding protein beta chain
<b>Expression Region</b>	16-258
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene encodes a member of a superfamily of proteins composed predominantly of tandemly arrayed short consensus repeats of approximately 60 amino acids. A single, unique beta-chain encoded by this gene assembles with seven identical alpha-chains into the predominant isoform of C4b-binding protein, a multimeric protein that controls activation of the complement cascade through the classical pathway. C4b-binding protein has a regulatory role in the coagulation system also, mediated through the beta-chain binding of protein S, a vitamin K-dependent protein that serves as a cofactor of activated protein C. The genes encoding both alpha and beta chains are located adjacent to each other on human chromosome 1 in the regulator of complement activation gene cluster. Alternative splicing gives rise to multiple transcript variants.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.