



Recombinant Human Protein FAM205B (FAM205B)

Product Code	CSB-MP720583HU
Abbreviation	FAM205B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63HN1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MHPMALHMVL PAHLPFLSPE VLRLLEVHVK KWMHFQRWGL PRRVEESLRQ LMPNPPLYQ PGNDQPVSFN LNNTSQVSLH RSETISLQTW CSCVAGQPIQ TFWVSEWSTM NPEQRHHCQQ TPNPMALALP SPALKALSGP HPQSGGQDND SGSDLQKYS QLFCGLPSLH SESLVATFMG SQGLPKIENV PKPPLKDPFL FNDLSFPQLL PKTSPQSAPP SSPLSPNWVS PSDHQRAQIN VPFLTAEYE ALEWHLLQRQ LQLQWGWPAALQRSQHTQCL MQHEPCGKAQ SPETTTASQT GKSISVLTRE LFFPEHARK LLEFHIQKQS IRHRWGLPQK IQQSIQLLLT STDQQTVSSS STALANVSIP QPVALEANGA CDVLSPIAAP VSIPRPHLLT QVKAILQSHI DSKCGQIHQG KIPACVHRSW DCRISGVLAV APFPCIPESQ FLVLQTASDP DLHHKVM PWM PTALDQQQQA LPGTVTEHPK LLRVLSVEAI EKLETTLRHK HLAFLSGLPA LYYVALPRAL APAVTSQSVI TEMEPSPVEI PAEPLI
Source	Mammalian cell
Target Names	FAM205BP
Protein Names	Recommended name: Protein FAM205B
Expression Region	1-556
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.