



Recombinant Burkholderia pseudomallei 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase (ispD)

Product Code	CSB-MP720612BXW
Abbreviation	ispD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q63T70
Product Type	Recombinant Protein
Immunogen Species	Burkholderia pseudomallei (strain K96243)
Purity	≥85% (SDS-PAGE)
Sequence	MTSRLFALIP CAGTGSRSRGS ALPKQYRTLA GRALLHYTLA AFDACSEFAQ TLVVISPDDA HFDARRFAGL RFAVRRCGGA SRQASVMNGL IQLAEFGATD ADWVLVHDAA RPGITPALIR TLIGALKDDP VGGIVALPVA DTLKRVPAGG DAIERTESRN GLWQAQTPQM FRIGMLRDAI RRAQLDGHDL TDEASAIWEA GHTPRVVQGS LRNFKVTYPE DFDLAEAILA QPARAS
Source	Mammalian cell
Target Names	ispD
Protein Names	Recommended name: 2-C-methyl-D-erythritol 4-phosphate cytidyltransferase EC= 2.7.7.60 Alternative name(s): 4-diphosphocytidyl-2C-methyl-D-erythritol synthase MEP cytidyltransferase Short name= MCT
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.