



Recombinant Mouse Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C (Pde1c)

Product Code	CSB-EP720665MO-B
Abbreviation	Pde1c
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64338
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MESPTKEIEE FESNSLKHLQ PEQIEKIWLR LRGLRKYKKT SQRLRSLVKQ LERGEASVVD LKKNLEYAAT VLESVYIDET RLLDTEDEL SDIQSDAVPS EVRDWLASTF TRQMGMMLRR SDEKPRFKSI VHAVQAGIFV ERMYYRRTSNM VGLSYPPAVI DALKDVDTWS FDFVSLNEAS GDHALKFIFY ELLTRYDLIS RFKIPISALV SFVEALEVGY SKHKNPYHNL MHAADVTQTV HYLKYTGVA NWLTELEIFA IIFSAAIHDY EHTGTTNMFH IQTRSDPAIL YNDRSVLENH HLSAAYRLLQ EDEEMNILVN LSKDDWREFR TLVIEMVMAT DMSCHFQQIK AMKTALQQPE AIEKPKALSL MLHTADISHP AKAWDLHHRW TMSLLEEFFR QGDREAELGL PFSPLCDRKS TMVAQSQVGF IDFIVEPTFT VLTDMTKIV SPLIDESSQT GGTGQRRSSL NSINSSDAKR SGVKSSGSDG SAPINNSVIP VDYKSFKATW TEVVQINRER WRAKVPKEEK AKKEAEEKAR LAAEEKQKEM EAKSQAEQGT TSKGEKKTSG EAKSQVNGTR KGDNPRGKNS KGEKAGEKQQ NGDLKDGNK ADKKDHSNTG NESKKTDTGK KRSHGSPAPS TSSTSRITLP VIKPPLRHFV RPAYASSSYA PSVPKKTDDH PVRKMLDQR IKMKKIQNIS HHWNKK
Source	E.coli
Target Names	Pde1c
Protein Names	Recommended name: Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C Short name= Cam-PDE 1C EC= 3.1.4.17
Expression Region	1-706
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.