



Recombinant Human Protein FRG2 (FRG2)

Product Code	CSB-YP720715HU
Abbreviation	FRG2
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q64ET8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGKGNEDSDL HCSSIQCSTD QPPFQQISFT EKGSDEKKPF KEKGKTAFSH SSEKHIQRQG SEPNNKENS EETKLKAGNS TAGSEPESSS YRENCRKRKM SSKDSCQDTA GNCPEKECSL SLNKKSRSST AVHNSEIQET CDAHHRGHSR ACTGHSKRHR SRALGVQTPS IRKSLVTSVR AMSEAVYQDL AQVWAQQIHS PLTCEQLTLL TRLRGPLCAQ VQTLYSMATQ AAYVFPAESW LVPATLPGPG ESALDREAHP FPGQEITETV SGSDEAKL
Source	Yeast
Target Names	FRG2
Protein Names	Recommended name: Protein FRG2 Alternative name(s): FSHD region gene 2 protein
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>