



Recombinant *Drosophila teissieri* Kinetochore protein Spc25 (Spc25)

Product Code	CSB-MP720716DMO
Abbreviation	Spc25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64EW3
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila teissieri</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MAIIPESSY ERRVKALYEK QIRMEALEGK FIKKVYKFNS NLLDVKEAVL RHQRKVGKQLQ KVVMERREEL EKRVSFMEEL AQELEATKLR NLAMKEQIKQ RKMIARQRKN EIMERIQTLS KTTGTYVNQE ALPARVKGVT VLRGDKRDQL IPFDLNATDA EGLNSLCQHL ESLNVDVSQW QQLVSLAMDV AMEARAPTP PKEVANCKSI IEIDLTSPPTS HQA
Source	Mammalian cell
Target Names	Spc25
Protein Names	Recommended name: Kinetochore protein Spc25
Expression Region	1-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.