



Recombinant Human WD repeat-containing protein 25 (WDR25)

Product Code	CSB-MP720719HU
Abbreviation	WDR25
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q64LD2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTARTLSLMA SLVAYDDSDS EAETEHAGSF NATGQQKDTS GVARPPGQDF ASGTLDPKA GAQPTKHGSC EDPGGYRLPL AQLGRSDWGS CPSQRLQWPG KEPQVTFPIK EPSCSSLWTS HVPASHMPLA AARFKQVKLS RNFPKSSFHA QSESETVGKN GSSFQKKKCE DCVVPYTPRR LRQRQALSTE TGKGKDVPEQ GPPAGRAPAP LYVGPVSEF IQPYLNHYK ETTVPRKVLV HLRGHRGPVN TIQWCPVLSK SHMLLSTSMK KTFKVVNAVD SGHCLQTYSL HTEAVRAARW APCGRRILSG GFDFALHLTD LETGTQLFSG RSDFRITTLK FHPKDHNIFL CGGFSSEMKA WDIRTGKVMR SYKATIQQTL DILFLREGSE FLSSTDASTR DSADRTIIAW DFRTSAKISN QIFHERFTCP SLALHPREPV FLAQTNGNYL ALFSTVWPYR MSRRRRYEGH KVEGYSVGCE CSPGGDLLVT GSADGRVLMY SFRTASRACT LQGHTQACVG TTYHPVLPSV LATCSWGGDM KIIWH
Source	Mammalian cell
Target Names	WDR25
Protein Names	Recommended name: WD repeat-containing protein 25
Expression Region	1-544
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.