



Recombinant Orgyia pseudotsugata multicapsid polyhedrosis virus Late expression factor 3 (LEF-3)

Product Code	CSB-MP720776OEU
Abbreviation	LEF-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65365
Product Type	Recombinant Protein
Immunogen Species	Orgyia pseudotsugata multicapsid polyhedrosis virus (OpMNPV)
Purity	≥85% (SDS-PAGE)
Sequence	MMAAKREHAD CGAEPARKRV KENYKRVTKG LLSKMTISLE NHLYYTFTFR LLNENKTEAY YGNLQCFKDL VEQECYDVSL NFVKTKYNER IEINEYSKCD AAIDENVAVK LCLTRADFEN EEIVNVLAKL KCVFKRLGAN NYKMFVDINM QDAGGAVFVQ QVECFANLKV LASAAKAFVK SPDNFNSLMD FYYKNTNTLF YVHGVRCQHT SKGQNEFLNW TAGPSTSLET PSNTDNEDIY NLVHSHSTNN ISRANKHLKS MQLSLFKAQEQ KINDNGKHSF SVQFKTLDSM DDDDDTKWHK CVYYVDSNGN KEDPNDTNAV QKLAMDFDQL ATCLADKLTAT IFVTADNA DASTMNLGL LKHDDEECEY QFL
Source	Mammalian cell
Target Names	LEF-3
Protein Names	Recommended name: Late expression factor 3
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.