



# Recombinant *Bacillus licheniformis* Energy-coupling factor transporter ATP-binding protein EcfA 1 (ecfA1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP720849BQU-B                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Abbreviation</b>      | ecfA1                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                          |
| <b>Uniprot No.</b>       | Q65P77                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Immunogen Species</b> | <i>Bacillus licheniformis</i> (strain ATCC 14580 / DSM 13 / JCM 2505 / NBRC 12200 / NCIMB 9375 / NRRL NRS-1264 / Gibson 46)                                                                                                                                                                                                                                                                                 |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Sequence</b>          | MNQERLIAAE DITFRYQKDA DRPALDHVSF HVSKGEWLAI VGHNGSGKST<br>LARALNGLIL PDEGSIKVGK IKLTEETVWD VRKKVGMVFQ NPDNQFVGTT<br>VRDDVAFGLE NSGVPREDMV ERVDWAVQQV NMQEFLDQEP<br>HHLSGGQKQR VAIAGVLAAR PDIMILDEAT SMLDPMGREE VLDTVRRLKD<br>QGMVTVISIT HDLNEAAKAD RIIVMNSGRK FAEGTPQEVF RLNQDLIKIG<br>LDLPFSFQLT QMLKESGLAL KGDHLTQEGL VNELWTLKSK M                                                                        |
| <b>Source</b>            | <i>E.coli</i>                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Target Names</b>      | ecfA1                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Energy-coupling factor transporter ATP-binding protein EcfA 1<br>Short name= ECF transporter A component EcfA 1 EC= 3.6.3.-                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-281                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                         |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                               |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                           |



of lyophilized form is 12 months at -20°C/-80°C.