



Recombinant Cytidylate kinase 2 (cmk2)

Product Code	CSB-MP720919BUI
Abbreviation	cmk2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q65ZV5
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MKIALSGKSG CGNTTVSSMI AKHYGLEFIN YTFHDIARER NIPFSEFYEK EIIGRDDYYW DMYLDKRLSL LSKKNNTVLA SRLAIWISKS ADLKIYLYAK MEVRAERIMT REGGMYSVDL SSTFNRDEND KKRYLSIYNI DIDDYSSKTD LVIDVTNINS NEVFELIRDE IDKRNLKN
Source	Mammalian cell
Target Names	cmk2
Protein Names	Recommended name: Cytidylate kinase 2 Short name= CK 2 EC= 2.7.4.25 Alternative name(s): Cytidine monophosphate kinase 2 Short name= CMP kinase 2
Expression Region	1-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.