



Recombinant GTPase Era (era)

| | |
|--------------------------|---|
| Product Code | CSB-MP720926BUI |
| Abbreviation | era |
| Storage | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |
| Uniprot No. | Q660L0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MKSGFAAILG RPSTGKSTLL NSICGHKISI ISPIPQTTRN KIKGIFTDDR GQIIFIDTPG FHLSKKKFNI AMMKNIHSSI GEVELILYII DIQDKPGEEE NKILEIKNK KIKFLVLLNK VDLKNTKIKE ITQFLKEKEI EDNNIIKISA EKKINTEELK NKIYENFSEG PLYYPQEYYT DQEINFRISE IIREKAIENL KEELPYSLYV DIDILENKKG SLFIRANIFV ANESQKGIIV GKNGKEIKSI GERARITIAK IFETKCNLFL QVKLKKNWNK EDKLIKRLIN |
| Source | Mammalian cell |
| Target Names | era |
| Protein Names | Recommended name: GTPase Era |
| Expression Region | 1-290 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |