



# Recombinant Serine--tRNA ligase (serS)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP720945BUI                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Abbreviation</b>      | serS                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                         |
| <b>Uniprot No.</b>       | Q662D8                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Immunogen Species</b> | Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBI) (Borrelia baviensis)                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Sequence</b>          | MLDLKFIRDN LDLVKKSIA RGLVLDIDKL VYLDDKRKKI ITQIGELNAK<br>RNENSSKMRE NLDNVFKSSL METGKILKKQ LIDLEEELKK ISFDFDLENK<br>RVPNISPVDV PIGNSEEDNF EIKKVGVPK FDFKPKDHLE LGRDLDLLDF<br>DRAREVSGSK FYYLKNEAVF LEIALINFSL NKLREKGFVD FITPDVAREF<br>IVDGIGFNPR GDESNIYKIE STDKYLVGTS EITLGGYYYN KIIDLTLPIK<br>MAGFSHCFRK EAGAYGQLSK GLYRVHQFSK VEMFCFCKAE ESQVIHDEFL<br>SIQEIQIFTEL EIPYRVNLIC SFDLGSPAYK KYDIEAWMPG RDGKGSGYGEV<br>TSTSNCTDYQ SRRLKIRYKD QDGQNKFAHM VNGTAIATTR VIISILENFQ<br>DEKGGVRIPK SLVKYTGFY IPFKN |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Target Names</b>      | serS                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Protein Names</b>     | Recommended name: Serine--tRNA ligase EC= 6.1.1.11 Alternative name(s): Seryl-tRNA synthetase Short name= SerRS Seryl-tRNA(Ser/Sec) synthetase                                                                                                                                                                                                                                                                                                                                                             |
| <b>Expression Region</b> | 1-425                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                         |