



Recombinant Peptide chain release factor 1 (prfA)

Product Code	CSB-EP720946BUI-B
Abbreviation	prfA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q662G9
Product Type	Recombinant Protein
Immunogen Species	Borrelia baviensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia baviensis)
Purity	≥85% (SDS-PAGE)
Sequence	MFLEKLSAT SKIKLIEEKL QDINLIKNQK KYSKIIKEYT YLEKINTKKI EYEKILSQIN DNKTILEKED QQEMKELIKQ ELIDLDKKKE DLEHQIKILL LPQDENDSKN IIIIEIRAGTG GEEAALFANN LYSMYIKYSE RKKWKTEIIN FNETELGGFK EIVFEIKGKD VFKKLKYESG VHRVQRIPIT ESNGLRQTSA ATVAVLPNIE ETEIDINEKD LRIDVYRSSG AGGQHVNTTD SAVRITHLPT GIVVQCQNER SQHKNKDQAM KILRARLYEF EDSKKQEQRS SNRKQQVGS DRSERIRTYN FPQNRITDHR ANITLYKLEE FMQGELDPLL DPLTIELQEQ KFKSNNI
Source	E.coli
Target Names	prfA
Protein Names	Recommended name: Peptide chain release factor 1 Short name= RF-1
Expression Region	1-357
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>