



Recombinant tRNA modification GTPase MnmE (mnmE)

Product Code	CSB-MP720947BUI
Abbreviation	mnmE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q662I5
Product Type	Recombinant Protein
Immunogen Species	Borrelia bavariensis (strain ATCC BAA-2496 / DSM 23469 / PBi) (Borrelia bavariensis)
Purity	≥85% (SDS-PAGE)
Sequence	MSNLFQRDDD IVALATPFLS SALCVIRSSG ASSIPKFSKI FSNHSALNSA PGNTIHYGYI LDNENCKVD EVVCLYRAP KSFTGQDSIE VIAHGSVIGI KKIIDLFLKS GFRMAEPGEF TLRSLAKKI ALTKAEALHE IIFATTNKTY SLAVTPLSGA LFANLDAIQR CLLNILSAVS AYLDYEVDDH EISIPFDLIL NSKAELKKLI NSYKVYEKID HGVTLVLAGS VNAGKSSLFN MFLKKDRSIV SSYPGTTDRDY IEATFELDGI LFNLFDTAGL KDADNFVERL GIEKSNSLIK EASLVIYVID VSSNLPREDF LFIDSNKSNS KILFVLNKID LNINKSTEEF VRSSVLNSSN LIMISTKNLE GIDILYDKIK TLISYERVEI GLDDIISSS RQMQLLEKAY ALVLDLLNKI DRQVSYDMLA FDAYEIINCL GEITGEVSSE DVLDNMFKNF CLGK
Source	Mammalian cell
Target Names	mnmE
Protein Names	Recommended name: tRNA modification GTPase MnmE EC= 3.6.-.-
Expression Region	1-464
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.