



Recombinant GTPase Der (der)

Product Code	CSB-MP720992YAH
Abbreviation	der
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q668A3
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MIPVIALVGR PNVGKSTLFN RLTHTRDALV ADFPGLTRDR KYGRAEVEGH EFIVVDTGGI DGTEDGVETK MAGQSLLAIE EADIVLFMVD ARAGLMPADQ GIAQHLSRE KATFLVANKT DGIDPDTATA DFYSLGLGEV HAI AASHGRG VTQLIEDVMA PYMDAEEPEV ELTDEEENAA YWAEQEAQGE DVPPEDPEDD FDPRTLPIKL AIVGRPNVGK STLTNRILGE DRVVVYDMPG TTRDSIYIPM TRDDREYILI DTAGVRKRGK ITETVEKFSV IKTLQAIEDS NVVLLVIDAR DGISDQDLSL LGFILNSGRS LVI AVNKWDG MTEEARAQVK DMLDLRLGFV DFARIHFISA LHGSGVGNLF ESVQEAYDCS TKRVGTSLLT RIMQMAEEDH QPPLVRGRRV KLKYAHAGGY NPPIVVIHGN QVTDLSDSYK RYLMNYFRRS LKVMGTPIRI QFKEGENPFA GKRNPLTPNQ MRKRKRLMSH LKKGK
Source	Mammalian cell
Target Names	der
Protein Names	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.