



Recombinant Tryptophan synthase beta chain (trpB)

Product Code	CSB-MP721009YAH
Abbreviation	trpB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66AK4
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTTLNPYFGE FGGMYVPQIL VPALKQLEDA FVSAQLDPEF QAAFQDLLKN YAGRPTALT CQNLTGKT KT KLYLKRELL HGGAHKTNQV LGQALLAKRM GKTEIIAETG AGQHGVASAL ACALLGLKCR IYMGAKDIER QSPNVFRMRL MGAEVIPVHS GSSTLKDACN EALRDWSGTY ETAHYMLGTA AGPHPYPTIV REFQRMIGEE TKAQILEKEG RLPDAVLACV GGGSNAIGMF ADFIDEPDVG LIGVEPAGLG IETGQHGAPL KHGKVGIYFG MKSPMMQTSD GQIEESYSIS AGLDFPSVGP QHAYLNSIGR ADYVSITDDE ALDAFKTLSC KEGIIPALES SHALAHALKM IKADPDKEQI LVVNLSGRGD KDIFTVHDIL KARGEI
Source	Mammalian cell
Target Names	trpB
Protein Names	Recommended name: Tryptophan synthase beta chain EC= 4.2.1.20
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.