



Recombinant tRNA (uracil (54)-C (5))-methyltransferase (trmA)

Product Code	CSB-EP721055YAH
Abbreviation	trmA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66G58
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MTPNILPIES YDHQLAEKSA RLKAMMLPFQ APEPEIFRSP ADHYRMRAEF RVWHDEDDLY HIMFDQQTQK RIRVEQFPVA SRLINRLMDA LMTAIRAEPL LRRKLFQIDY LSTLSGKLIA SLLYHRQLDE EWQQKALELR DQLRAQGFDL QLIGRAAKTK IMLDHDYIDE VLPVAGREMI YRQVENSFTQ PNAAVNIHML EWALDVTQGA TGDLELYCG NGNFSALAR NFERVLATEI AKPSVAAAQY NIAANNIDNV QIIRMSAEFF TQAMQGVREF NRLKGIDLGS YNCETIFVDP PRSGLDHETV KLVQAYPRIL YISCNPETLC ANLEQLQHTH KISRLALFDQ FPYTHHMECG VLLEKRH
Source	E.coli
Target Names	trmA
Protein Names	Recommended name: tRNA (uracil(54)-C(5))-methyltransferase EC= 2.1.1.35 Alternative name(s): tRNA(m5U54)-methyltransferase Short name= RUMT
Expression Region	1-367
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.