



Recombinant Triosephosphate isomerase (tpiA)

Product Code	CSB-EP721057YAH
Abbreviation	tpiA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66GA1
Product Type	Recombinant Protein
Immunogen Species	Yersinia pseudotuberculosis serotype I (strain IP32953)
Purity	≥85% (SDS-PAGE)
Sequence	MRHPLVMGNW KLNGSTHVMN ELIAGLRKEL STVDGCGVAI APPAIYLNQA KHELGSRIA LGAQNV DVNL SGAFTGETSA EMLKDIGAQY IIIGHSERRT YHQESDELIA KKFGVLKEIG LIPVLCIGES EAENEAGQTE AVCAKQLDAV LNTLGVKAFE GAVIAYEPIW AIGTGKSATP AQAQAVHKFI RDHIAKQDAA VAAQVIIQYG GSVNDKNAAE LFTQPDIDGA LVGGASLKAD AFAVIVKAAA KAKKA
Source	E.coli
Target Names	tpiA
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.