



Recombinant *Macaca fascicularis* Protein FAM170A (FAM170A)

Product Code	CSB-EP721100MOV
Abbreviation	FAM170A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q66LM5
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	≥85% (SDS-PAGE)
Sequence	MKRRQKRKHL ENEESQETAE KGGGIQRIHR DSPQPQSPLA QVQERGETPP RSQHVLSLFH SSYKTCVSSL CVNKEERGMK IYYMQVQMNK GVAVSWETEE TLESLEKQPR MEEVTLSEVV RVGTPPSDVS TRNLLSDSEP SGEEKEHEER TESDSLPGSP TVEDTPRAKT PDWLVTMENG FRCMACCRVF TTMEALQEHV QFGIREGFSC HVFHLTMAQL TGNMESESTQ DEQEEENGNE KEEEEKPEAK EEEGQPTEED LGLRRSWSQC PGCVFHSPKD RNS
Source	<i>E.coli</i>
Target Names	FAM170A
Protein Names	Recommended name: Protein FAM170A Alternative name(s): Zinc finger domain-containing protein Zinc finger protein ZNFD
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.