



Recombinant Wheat dwarf virus Replication-associated protein (C1/C2)

Product Code	CSB-MP721122WAH
Abbreviation	C1/C2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67622
Product Type	Recombinant Protein
Immunogen Species	Wheat dwarf virus (isolate Sweden) (WDV)
Purity	≥85% (SDS-PAGE)
Sequence	MASSAPRFR VYSKYLFLTY PECTLEPQYA LDSLRLLNK YEPLYIAAVR ELHEDGSPHL HVLVQNKLRA SITNPNALNL RMDTSPFSIF HPNIQAAKDC NQVRDYITKE VDSDVNTAEW GTFVAVSTPG RKDRDADMKQ IIESSSSREE FLSMVCNRFP FEWSIRLKDF EYTARHLFPD PVATYTPEFP TESLICHETI ESWKNEHLYS ESPGRHKSIY ICGPTRTGKT SWARSLGTHN YYNSLVDFTT YDVNAKYNII DDIPFKFTPN WKCFVGAQRD FTVNPKYGKR KVIRGGIPCI ILVNPDEDWL KDMTPEQSDY MYSNTVVHYM YEGETFINYS FASGEDVTAS Q
Source	Mammalian cell
Target Names	C1/C2
Protein Names	Recommended name: Replication-associated protein Short name= Rep EC= 2.7.7.- EC= 3.1.21.-
Expression Region	1-351
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.