



Recombinant Bifunctional protein GlmU (glmU)

Product Code	CSB-MP721145SRW
Abbreviation	glmU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67JC8
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MSDITAVLLA AGHGTRMKSD LIKVMHPLAG KPMIGHIVDN VRRAGLEDIV VVVGYQQERI REYLGDRVRY AVQSEQLGTG HAVLQAAGLI DETEGGHVLV MYGDNPFIGP ELIERLIRAH VEADAAASLL TAEADPGAL GRILRDPATG AFLGSVEYKD ATPEQRRIRE IWTGVAVFRR AGFTALLNRL DRNNAQGEYY LPQVWEILLQ RGEKVQALLL ASEEDALAPN DRVELARAEA RLRRQINERH MRNGVTIINP DATYIDEDVE IGRDTVWPF TFIHGKTVIG PHCKIGPMTT IVSSTVAEGC VVEQSVVEES YVGPGRIGP MAHLRPGCEL EGAAEIGNYA ELKKAKVGRG VKCHHHSYLG DATIGAGANI GAGTITANYN GVEKFRTEIG SGAFIGNVN LIAPITVGDG ALIAAGSTVG PRLEIPADAL VVERAQAVIK EGRAASLKEA WRQRKMNREG T
Source	Mammalian cell
Target Names	glmU
Protein Names	Recommended name: Bifunctional protein GlmU Including the following 2 domains: UDP-N-acetylglucosamine pyrophosphorylase EC= 2.7.7.23 Alternative name(s): N-acetylglucosamine-1-phosphate uridyltransferase Glucosamine-1-phosphate
Expression Region	1-471
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.