



Recombinant 3-isopropylmalate dehydrogenase (leuB)

Product Code	CSB-MP721153SRW
Abbreviation	leuB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q67JY2
Product Type	Recombinant Protein
Immunogen Species	Symbiobacterium thermophilum (strain T / IAM 14863)
Purity	≥85% (SDS-PAGE)
Sequence	MAQPTYHIAC LPGDGIGPEV TRGAVTVLQA AAAAYGFRLE FSEYLVGGAA YDAVGTPFPD ETRDACDRAD AILFGAVGGP RYEGLPWDLR PEAGLLAIRK RYDLYANLRP ILLYPPLKDA SPLKNEIIGG GVDFIIVREL VGGIYFGEPR GIETLPDGTR RGVNTEVYTD AEIARIARMG FEIARGRRRR LTSVDKGNVM EAGKLWRTVV DEVAREFPDV EVEHLLADNA AMQILRRPGD FDVMLASNLF GDFLSDEAAM LTGSIGLLPS ASLGAARNRF GLPKGFYEPI HGSAPDIAGQ DRANPLAAIL CGAMLLRHSL GREDAARAVE QAVAAVLEEG LRTADITVPG TAVLGTAAMA RAVADRVH
Source	Mammalian cell
Target Names	leuB
Protein Names	Recommended name: 3-isopropylmalate dehydrogenase EC= 1.1.1.85 Alternative name(s): 3-IPM-DH Beta-IPM dehydrogenase Short name= IMDH
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.