



# Recombinant GTPase Der (der)

<b>Product Code</b>	CSB-BP721167SRW
<b>Abbreviation</b>	der
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q67NS9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Symbiobacterium thermophilum (strain T / IAM 14863)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MKPIVAIVGR PNVGKSTLFN RLTQSRHAIV EDQPGVTRDR LYADTEWNGR TLTLVDTGGI QLDKEGDTIE AHVTRQAELA IREADVIIFV VDVTGVTAP DLEVADLLRR QRKPVIVAVN KVENLKREDE ALEFWALGLE PLINVSAEHG LGTGDLLDAV VAALPDLSEP EPEEGPVRV AVIGRPNVGK SSLVNAILGE ERVIVSDVPG TTRDAIDVLV ERGEDKFLLI DTAGMRRKAR VEEAVERYSV MRALRAVERA QVVLIVIDAQ DGVTEQDQRI AGYAHENGKA CIVVVKWDL IEKDDRTMAK MTEEVRMRLA FMDYAMIHFV SAKTRARVHR LLPLIKEAAA NHARRISTRE LNDLVREAVA LNPPSPDKGR RLKIFYATQP HVSPPGFVFF VNDSELVHFS YQRYLENQLR QTYAFEGTPI NLYFRTREKA KLGERPVRVR RVLAGGQTRE IRRAARKRST E
<b>Source</b>	Baculovirus
<b>Target Names</b>	der
<b>Protein Names</b>	Recommended name: GTPase Der Alternative name(s): GTP-binding protein EngA
<b>Expression Region</b>	1-471
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.