



# Recombinant Arabidopsis thaliana Beta-glucosidase 8 (BGLU8)

<b>Product Code</b>	CSB-BP721193DOA
<b>Abbreviation</b>	BGLU8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q67XN2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	FTRNDFPE DFLFGAGTSA YQWEGAANED GRTPSVWDTT SHCYNGSNGD IACDGYHKKYK EDVKLMAEMG LESFRFSISW SRLIPNGRGR INPKGLLFYK NLIKELRSHG IEPHVTLYHY DLPQSLEDEY GGWINHKIIE DFTAFADVCF REFGEDVKLW TTINEATIFA FAFYGKDVMY GNCTTGNYCM ETYIAGHNML LAHASASNLY KLKYKSKQRG SIGLSIFALG LTPYTNSKDD EIATQRAKAF LYGWMLKPLV FGDYPDEMKR TLGSRLPVFS EESEQVKGS SDFVGIIHYT TVYVTNQPAP YIFPSSTNKD FFTDMGAYII STGNSSSFVF DAVPWGLEGV LQHIKHRYNN PPIYILENGS PMKHDSMLQD TPRVEYIQAY IGAVLNAIKS GSDTRGYFVW SLIDLFEVQV GYKSSFGMYV VNFSDPGRKR SPKLSASWYT GFLNGTIDVA SQDMTQLQRN FSGSSSL
<b>Source</b>	Baculovirus
<b>Target Names</b>	BGLU8
<b>Protein Names</b>	Recommended name: Beta-glucosidase 8 Short name= AtBGLU8 EC= 3.2.1.21
<b>Expression Region</b>	23-497
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.