



Recombinant *Oryza sativa* subsp. japonica Protein DEHYDRATION-INDUCED 19 (DI19-1)

Product Code	CSB-EP721202OFG
Abbreviation	DI19-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q688X9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSEHWISRL AAAKRFYAAQ LGHADRAGME EVDMDEEVRP EFACPYCYED HDVVSLCAHL EEEHPFEPHA APCPICSDKI AKDMLNHITV QHGYLFKNRR RLRRFVIPGS QALSLLSRDL REAHLQVLLG GGGHRSNNSS NTTNISADPL LSSFGLSFPT SDTEETSKPP ISIPDDASVI KETPAQPWDS SIDSSLTREE REQKRKQASV RATFVQDLLL TTLFGD
Source	E.coli
Target Names	DI19-1
Protein Names	Recommended name: Protein DEHYDRATION-INDUCED 19 Alternative name(s): OsDi19
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.