



Recombinant Rickettsia typhi Tryptophan--tRNA ligase (trpS)

Product Code	CSB-BP721266RNE
Abbreviation	trpS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q68WR2
Product Type	Recombinant Protein
Immunogen Species	Rickettsia typhi (strain ATCC VR-144 / Wilmington)
Purity	≥85% (SDS-PAGE)
Sequence	MKKTVLSGVQ TTGALHLGNY LGSIKNWVKM QEEYNCFFFL ADLHAITIDI KTSELNNSIM EVLAIYLAAG LNPDKVTIFA QSMVKEHAEL SWLLNCVTPL GWLKRMTQFK DKAGSDQCKA CLGLFAYPIL MAADILYKA DIVPVGEDQK QHLELTRDIA EVINRRFDKE ILKVPDILIS ETGTRIMSLR NGLKKMSKSD IADFSRINLK DSNDLIYQKI KKAKTDSLFS ISYDKKTRPE ISNLLDIYKS LSKESMEKII NNYQNQGFAK FKEDLAEIII TNLQPIRDKC LELMNDKEYL LKILHKGAEK AKIRASETVN EVKKQFGFII
Source	Baculovirus
Target Names	trpS
Protein Names	Recommended name: Tryptophan--tRNA ligase EC= 6.1.1.2 Alternative name(s): Tryptophanyl-tRNA synthetase Short name= TrpRS
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.