



Recombinant *Oryza sativa* subsp. japonica Probable pyridoxal biosynthesis protein PDX1.1 (PDX11)

Product Code	CSB-BP721311OFG
Abbreviation	PDX11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q69LA6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MATDGTGVVT VYGSNGTNGAA LLEPSNHKSA TFSVKVGLAQ MLRGGVIMDV VTPEQARIAE EAGACAVMAL ERVPADIRAQ GGVARMSDPG LIRDIKRAVT IPVMAKARIG HFVEAQILEA IGVDYVDESE VLT LADDAHH INKHNFRVFP VCGCRDLGEA LRRIREGAAM IRTKGEAGTG NVVEAVRHVR SVMGDIRALR NMDDDEVFSY AKRIAAPYDL VMQTKQLGRL PVVQFAAGGV ATPADAALMM QLGCDGVFVG SGIFKSGDPA RRARAIVQAV THYSDPKILA EVSSGLGEAM VGINLSDPKV ERFAARSE
Source	Baculovirus
Target Names	PDX11
Protein Names	Recommended name: Probable pyridoxal biosynthesis protein PDX1.1
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.